TRANSMIT READING TO YOUR MOBILE DEVICE

This monitor syncs your blood pressure readings to the **Homedics Health+** app that is free to download on the App Store® and on Google Play™. Make sure you have the **Homedics Health+** app downloaded, and open on your mobile device before trying to transmit your blood pressure measurements. Before attempting to sync the blood pressure monitor with your mobile device, make sure Bluetooth® is turned ON on both your mobile device and the monitor.



- Transmit Reading**  
There are 2 ways to activate the Bluetooth® function. **Automatically activate:** When measurement completed, the device activates the Bluetooth® function automatically, and Bluetooth® symbol  will flash on the screen. **NOTE:** The app must be open for the data to transmit successfully. **Manually activate:** Under standby mode, user can press the  button for 3 seconds to activate the Bluetooth® function, and the Bluetooth symbol  will flash on the screen. **NOTE:** The app must be open for the data to transmit successfully.
- If the monitor is connected successfully to your smart device, the Bluetooth® symbol  will appear on the screen.
  - If the monitor cannot be connected to paired Bluetooth® - enabled device over 45 seconds, LCD will display Error message "E4" and Bluetooth® will be turned off.

- NOTE:**
- The measurement is stored in the monitor's memory even if the reading is not transmitted to your mobile device.
  - The **Homedics Health+** app must be open on your mobile device in order to transmit your measurements.
  - Only new readings will be accepted by the App.
  - This monitor can only pair up with one Bluetooth-enabled device at a time.
  - To ensure readings transmitted to the App have the correct date and time, it is important that the correct date and time is set on the blood pressure monitor before taking measurements.
  - Measurements transmitted to the App cannot be edited.
  - Bluetooth® compatibility with blood pressure monitor for Bluetooth-enabled device is:
    - Bluetooth® 4.2 for Android 6.0 or above
    - Bluetooth® 4.2 for iOS 7.0 or above

RECALLING VALUES FROM MEMORY

This monitor can be used by two individuals. Each user can store up to 120 measurements. This monitor features a memory mode to provide you with a variety of options for reviewing your measurement history, including AM/PM averaging. Reviewing your morning (AM) and nighttime (PM) measurements can provide important information about how your blood pressure changes throughout the day.

1. Press the **USER-SELECT**  button to choose User 1 or User 2.
  2. Press the  button to access the memory.
  3. The monitor will first display the calculated average applied to the last 3 measurements.
- NOTE:**
- If there is no data stored, nothing (except month, date, and time) will appear on the display.
  - If the records are less than 3, it will display the latest record.
4. Press the  button again and the monitor will display the average of the last 3 AM readings.
  5. Press the  button again to display the average of the last 3 PM readings.
  6. Every new press of the  button will recall a previous reading. The latest reading will be recalled first.
  7. To exit memory mode, press **START/STOP** button, and switch to sleep mode.

| Symbol                                                                             | Status                                                         |
|------------------------------------------------------------------------------------|----------------------------------------------------------------|
|  | An average of last 3 measurements.                             |
|  | An average of last 3 morning measurements. (4:00AM - 12:00PM)  |
|  | An average of last 3 nighttime measurements. (6:00PM - 2:00AM) |

**NOTE:** If the wrist positioning guide is ON before the measurement and the monitor is not in the ideal position,the  icon will be displayed with the saved readings.

CLEARING VALUES FROM MEMORY

1. Press the **USER-SELECT**  button to choose User 1 or User 2.
  2. Press the  button to access the memory recall mode.
  3. Press and hold the **DATE/TIME SET**  and  buttons at the same time and the data for the selected user will be erased automatically.
  4. To confirm deletion, press the  button and no data should appear.
- NOTE:** Once deleted, your readings cannot be restored.

IMPORTANT NOTES REGARDING YOUR BLOOD PRESSURE MEASUREMENT

- It is suggested that you take your measurements at the same time each day and use the same wrist for consistency.
- Users should wait a minimum of 5 minutes before taking additional measurements. More time may be necessary depending upon your physiology.
- The measurement results that users receive are for reference only. If users have any blood pressure concerns, please consult a physician.
- Once inflation reaches 300 mmHg, the unit will deflate automatically for safety reasons.
- This device is not suitable for people with arrhythmias.
- This device may have difficulty determining the proper blood pressure for pregnant women and for users with irregular heartbeat, diabetes, poor circulation of blood, kidney problems, or for users who have suffered a stroke.

CARE, MAINTENANCE, & CLEANING

- Clean the blood pressure monitor body and cuff carefully with a slightly damp, soft cloth. DO NOT press. DO NOT wash cuff or use chemical cleaner on it. Never use thinner, alcohol, or petrol (gasoline) as cleaner.
- Make sure the cuff is completely dry before using.
- Leaky batteries can damage the unit. Remove the batteries when the unit will not be used for a long time.
- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- If the unit is stored near freezing, allow it to acclimate to room temperature before use.
- This blood pressure monitor is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems with this device, please contact Homedics Consumer Relations (contact information on warranty page).
- DO NOT immerse the unit in water as this will result in damage to the unit.
- DO NOT subject the monitor or cuff to extreme temperatures, humidity, moisture, or direct sunlight. Protect from dust.
- DO NOT fold the cuff tightly.
- DO NOT disassemble the monitor or cuff. If in need of repair, refer to the warranty section of this manual.
- DO NOT subject the monitor to extreme shocks (DO NOT drop on floor).
- DO NOT inflate the cuff unless wrapped around wrist.
- DO NOT wrap the cuff around body parts other than your wrist.
- DO NOT drop or insert any object into any opening.
- This monitor may not meet its performance specifications if stored or used outside of these temperature and humidity ranges:

**Storage/Transportation Environment Temperature:** -13°F ~ 158°F(-25°C ~ 70°C)  
**Humidity:** Less than 93% RH  
**Operation Environment Temperature:** 41°F ~ 104°F(5°C ~ 40°C)  
**Humidity:** 15 ~ 93% RH  
**Atmospheric Pressure:** 700hPa-1060hPa

POTENTIAL FOR ELECTROMAGNETIC INTERFERENCE

To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipment, DO NOT use the device near a cell phone or microwave oven. For most wireless communication devices, it is recommended to maintain a distance of 10.8 feet (3.3m) in order to avoid electromagnetic interference.

**SUPPLIER'S DECLARATION OF CONFORMITY**  
**Product Description:** Blood Pressure Monitor  
**Model Number:** BPW-800BTA  
**Trade Name:** Homedics  
**FCC Compliance Statement**  
This device complies with Part 15 of the FCC Rules. This device contains license-exempt transmitter/receiver that complies with Innovation, Science and Economic Development Canada's license-exempt RSS(s). Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation.  
**U.S. Contact Information**  
**Company:** Homedics, LLC.  
**Address:** 3000 N. Pontiac Trail, Commerce Township, MI 48390  
8:30am-7:00pm EST Monday-Friday 1-800-466-3342  
Homedics is not responsible for any radio or TV interference caused by unauthorized modifications to this equipment. Such modifications could void the user authority to operate the equipment.

This equipment has been tested and found to comply with the Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

CAN ICES-003(B)/NMB-003(B)

RF Exposure Statement

This equipment complies with FCC and ISED RF radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operated in conjunction with any other antenna or transmitter.

This equipment essential performance per IEC 80601-2-30 additional essential performance requirements:

- 201.12.1.102 Limits of the error of the manometer from environmental conditions Over the temperature range of 5°C to 40°C (41°F ~ 104°F) and the relative humidity range of 15% to 93% (non-condensing), the maximum error for the measurement of the CUFF pressure at any point of the NOMINAL measurement range shall be less than or equal to ± 3 mmHg (± 0.4 kPa) or 2 % of the reading, whichever is greater.
- 201.12.1.107 Limits of the change in error of the blood pressure determination The laboratory limits of the change in error of the BLOOD PRESSURE DETERMINATION of the AUTOMATED SPHYGMOMANOMETER shall be less than 3 mmH (0.4 kPa).

ELECTROMAGNETIC COMPATIBILITY (EMC)

**Guidance and manufacturer's declaration – electromagnetic emissions**  
The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

| Emissions Test                                        | Compliance     | Electromagnetic Environment - Guidance                                                                                                                                                                                   |
|-------------------------------------------------------|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF emissions CISPR 11                                 | Group 1        | RF energy is used only to maintain device's operation. Therefore, its RF emissions are so low that it's not likely to cause any interference in nearby                                                                   |
| RF emissions CISPR 11                                 | Group B        |                                                                                                                                                                                                                          |
| Harmonic Emissions IEX 61000-3-2                      | Not Applicable | The device is suitable for use in all establishments, including domestic establishments, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Not Applicable |                                                                                                                                                                                                                          |

**Guidance and manufacturer's declaration – electromagnetic immunity**  
The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

| Immunity Test                                | IEC 60601 Test Level                              | Compliance Level              | Electromagnetic Environment– Guidance                                                                                                                              |
|----------------------------------------------|---------------------------------------------------|-------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD) IEC 61000-4-2  | ± 8 kV contact discharge<br>± 15 kV air discharge | ± 8 kV contact<br>± 15 kV air | In the case of air discharge testing, the climatic conditions shall be within the following ranges: Ambient Temperature: 15°C ~ 35°C Relative Humidity: 30% ~ 60%. |
| Power frequency magnetic field IEC 61000-4-8 | 30 A/m<br>50 or 60 Hz                             | 30A/m<br>50 or 60 Hz          | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                          |

**Guidance and manufacturer's declaration – electromagnetic immunity**  
The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

| Immunity Test              | IEC 60601 Test Level                                                                                                                                                                                                                                                                 | Compliance Level                                                                                                                                                                                                                                                                     | Electromagnetic Environment– Guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted RF IEC 61000-4-6 | 3V rms At 0.15-80 MHz 6V rms At ISM & RadioAmateur Freq                                                                                                                                                                                                                              | Not Applicable                                                                                                                                                                                                                                                                       | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Radiated RF IEC 61000-4-3  | 10 V/m at 80-2700 MHz<br>AM Modulation And 9-28V/m at 385-6000MHz- Pulse Mode and other Modulation. The system shall be tested as specified in IEC60601-1-2 Table 9 for proximity fields from RF wireless communications equipment using the test methods specified in IEC 61000-4-3 | 10 V/m at 80-2700 MHz<br>AM Modulation And 9-28V/m at 385-6000MHz, Pulse Mode and other Modulation. The system shall be tested as specified in IEC60601-1-2 Table 9 for proximity fields from RF wireless communications equipment using the test methods specified in IEC 61000-4-3 | Recommended separation distance<br>Considering to reduce the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:<br>$E=6/d\sqrt{P}$<br>where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVELS in V/m. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,a should be less than the compliance level in each frequency range." Interference may occur in the vicinity of equipment marked with the following symbol:  |

**NOTE 1:** At 80 MHz and 800 MHz, the higher frequency range applies.  
**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a: Field strengths from fixed trans-mitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast can-not be predicted theoretically with accuracy. To assess the electro-magnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.  
b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Test specifications for enclosure port immunity to RF wireless communications equipment.

| Test frequency (MHz) | Modulation                            | Immunity Test Level (V/m) |
|----------------------|---------------------------------------|---------------------------|
| 385                  | Pulse modulation 18 Hz                | 27                        |
| 450                  | FM± 5 kHz deviation 1kHz sine         | 28                        |
| 710                  |                                       |                           |
| 745                  | Pulse modulation 217 Hz <sup>a)</sup> | 9                         |
| 780                  |                                       |                           |







| Test frequency (MHz) | Modulation                            | Immunity Test Level (V/m) |
|----------------------|---------------------------------------|---------------------------|
| 810                  | Pulse modulation 18 Hz <sup>a)</sup>  | 28                        |
| 870                  |                                       |                           |
| 930                  |                                       |                           |
| 1720                 | Pulse modulation 217 Hz <sup>a)</sup> | 28                        |
| 1845                 |                                       |                           |
| 1970                 |                                       |                           |
| 2450                 | Pulse modulation 217 Hz <sup>a)</sup> | 28                        |
| 5240                 | Pulse modulation 217 Hz <sup>a)</sup> | 9                         |
| 5500                 |                                       |                           |
| 5785                 |                                       |                           |

**NOTE:** If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a). The carrier shall be modulated using a 50% duty cycle square wave signal. b). As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

TROUBLESHOOTING

If any abnormality arises during use, please check the following points:

| Symbols/ Symptoms                                                                                                                                                                                                     | Conditions/Causes                                                                                                                                                                             | Indication/Correction                                                                                                                                                            |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  Unit does not turn on when START/STOP button is pushed.                                                                           | Worn out batteries                                                                                                                                                                            | Replace them with 2 new AAA (LR03) alkaline batteries.                                                                                                                           |
|                                                                                                                                                                                                                       | Battery polarities have been positioned incorrectly.                                                                                                                                          | Reinsert the batteries in the correct positions.                                                                                                                                 |
|  Measuring Error Symbol appears when blood pressure value displayed is excessively low or high.                                    | Cuff has been placed incorrectly                                                                                                                                                              | Wrap the cuff properly so that it is positioned correctly.                                                                                                                       |
|  Measuring Error Symbol                                                                                                           | Air circuit abnormality. Cuff tube may not be plugged into monitor correctly.                                                                                                                 | Check cuff connection. Measure again.                                                                                                                                            |
|  Measuring Error Symbol                                                                                                          | Inflation pressure exceeding 300 mmHg.                                                                                                                                                        | Switch the unit off, then measure again.                                                                                                                                         |
|  Measuring Error Symbol                                                                                                          | Can't determine blood pressure measurement data.                                                                                                                                              | Wrap the cuff snugly so that it is positioned correctly. If you have any question about the cuff wearing and/or measurement result, please consult your healthcare professional. |
|                                                                                                                                                                                                                       | Cuff is worn improperly                                                                                                                                                                       |                                                                                                                                                                                  |
|  Measuring Error Symbol                                                                                                          | If the device cannot be connected to paired Bluetooth-enabled device over 45 seconds, LCD will display (Error message "E4" and Bluetooth® will be turned off.                                 | Please press  button for 3 seconds to start Bluetooth® function                             |
|                                                                                                                                                                                                                       |                                                                                                                                                                                               |                                                                                                                                                                                  |
|  Excessive Body Motion Detector<br><br>Notice: The measured blood pressure reading may not be accurate if the icon is displayed. | Body movement during measurement, especially, the movement on the wrist the blood pressure monitor is worn on. e.g. Talking, moving or shaking of the arm with the cuff on while measurement. | Measure again. Keep steady during measurement.                                                                                                                                   |
|                                                                                                                                                                                                                       | Cuff is worn improperly, or the shape of the wrist is unusual, excessive gap might be exist between the cuff and the wrist.                                                                   | Wrap the cuff properly and keep steady. Measure again. If you have any question about the cuff wearing and/or measurement result, please consult your healthcare professional.   |
|  Cuff Fitness detection Symbol                                                                                                   | The cuff was wrapped incorrectly (for example too loosely or too tightly).                                                                                                                    | Please reference "Applying the Cuff" section to wrap the cuff correctly.                                                                                                         |
|                                                                                                                                                                                                                       | BPM cannot communicate with Bluetooth-enabled device                                                                                                                                          | Please reference "Applying the Cuff" section to wrap the cuff correctly.                                                                                                         |
|                                                                                                                                                                                                                       | Bluetooth® function is not turned on.                                                                                                                                                         | See the "Bluetooth® Operation" section to turn on Bluetooth® function.                                                                                                           |
|                                                                                                                                                                                                                       | The distance between BPM and Bluetooth-enabled device is out of transmitting range.                                                                                                           | Please make sure the acceptable distance (<10 meters) with each other.                                                                                                           |
|                                                                                                                                                                                                                       | Use an incompatible Bluetooth-enabled device.                                                                                                                                                 | Please refer to "Bluetooth® Operation" and "RF Specification" sections.                                                                                                          |
|                                                                                                                                                                                                                       | Use non-Bluetooth-enabled device.                                                                                                                                                             |                                                                                                                                                                                  |
|                                                                                                                                  | Unexpected loss of electrical/mechanical integrity.                                                                                                                                           | Reinsert the batteries and try again.                                                                                                                                            |

If the unit still does not work, contact Homedics Consumer Relations. Under no circumstance should you disassemble or attempt to repair the unit by yourself. Contact information for Homedics Consumer Relations Department can be found in the Warranty page.

SPECIFICATIONS



|                                 |                                          |
|---------------------------------|------------------------------------------|
| Model Number                    | HL158B1-BD/ BPW-800BT-WT / BPW-800BTA-WT |
| Measurement Method              | Oscillometric                            |
| Rated Range of Cuff Pressure    | 0 ~ 300 mmHg                             |
| Rated Range of Determination    | 40 ~ 280 mmHg                            |
| Measurement Range of Heart Rate | Pulse: 40 ~ 199 Beats / Minute           |
| Accuracy                        | Pressure: ± 3 mmHg<br>Pulse: ± 5 % Max.  |
| Inflation                       | Automatic Inflation (Air Pump)           |

SPECIFICATIONS (CONT.)

|                                                    |                                                                                                                |
|----------------------------------------------------|----------------------------------------------------------------------------------------------------------------|
| Deflation                                          | Automatic (Passive Exhaust Valve)                                                                              |
| Display                                            | Liquid Crystal Display                                                                                         |
| Memory                                             | 240 Memory Total for 2 Users                                                                                   |
| Unit Dimensions                                    | 73.79 mm x 60.15 mm x 22.94 mm (LxWxH)<br>2.91 x 2.37 x 0.90 inch (LxWxH)                                      |
| Unit Weight                                        | 69 g ± 3 g (Without Cuff & Batteries)                                                                          |
| Cuff Size                                          | 5.3 ~ 7.7 inch (approx.135 ~ 195 mm)                                                                           |
| Storage/Transportation Environment                 | Temperature: -25°C ~ 70°C (-13°F ~ 158°F)<br>Humidity: ≤ 93% R.H.                                              |
| Operation Environment                              | Temperature: 5°C ~ 40°C (41°F ~ 104°F)<br>Humidity: 15% ~ 93% R.H.<br>Atmospheric pressure: 700 hPa ~ 1060 hPa |
| Power Supply                                       | DC 3V AAA "LR03" (1.5V)<br>Alkaline Battery x 2                                                                |
| Battery Life                                       | Approx. 250 Measurements                                                                                       |
| Shelf life (battery)                               | 3 years (Temperature: 20 ± 2°C;<br>Relative humidity: 65 ± 20% RH)                                             |
| Product Life                                       | 5 Years (4 times per day)                                                                                      |
| Sleeping Mode                                      | Without any operation for 1 minute, device automatically shuts off.                                            |
| Accessories                                        | Instruction Manual, 2 AAA (LR03) Alkaline Batteries, Storage Pouch                                             |
| RF Type                                            | Bluetooth® 5.0 BLE                                                                                             |
| System requirement of the Bluetooth-enabled device | Bluetooth® 4.2 for Android 6.0 or above<br>Bluetooth® 4.2 for iOS 7.0 or above                                 |

\* The contents of this manual and the specifications of the device covered by this manual are subject to change for improvement without notice.

EXPLANATION OF SYMBOLS

|                                                                                       |                                                                                                                                                                                                                                                                                                                                                                                                  |
|---------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | Follow instruction for use.                                                                                                                                                                                                                                                                                                                                                                      |
|    | TYPE BF Applied Part.                                                                                                                                                                                                                                                                                                                                                                            |
|   | To avoid inaccurate results caused by electromagnetic interference.<br><b>WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm (12 inches) to any part of the device. Otherwise, degradation of the performance of this equipment could result.</b>                                             |
|  | Waste of electrical and electronic equipment (VEEE). Discard the used product to the recycling collection point according to local regulations.                                                                                                                                                                                                                                                  |
|  | Serial number<br> YMMMXXXXX                                                                                                                                                                                                                                                                                 |
|  | Ingress Protection Rating.First characteristic numeral - Degree of protection against access to hazardous parts and against solid foreign objects N1=2 (Protected against solid foreign objects of 12.5 mm ø and greater). Second characteristic numeral Degree of protection against ingress of water. N2=2 (Protected against vertically falling water drops when ENCLOSURE tilted up to 15°). |
|  | Non-ionizing electromagnetic radiation.                                                                                                                                                                                                                                                                                                                                                          |

Homedics

For service or repair, DO NOT return this unit to the retailer. Contact Homedics Consumer Relations at:  
  
Email: cservice@homedics.com  
  
Phone: 1-800-466-3342  
Business Hours: 8:30am-7pm EST Monday-Friday  
  
Distributed by: Homedics, LLC  
3000 N Pontiac Trail  
Commerce Township, MI 48390  
Printed in China or Vietnam

**LIMITED 5-YEAR WARRANTY**  
Homedics sells its products with the intent that they are free of defects in manufacture and workmanship for a period of 5 years from the date of original purchase, except as noted below. Homedics warrants that its products will be free of defects in material and workmanship under normal use and service. This blood pressure monitor meets the simulated measurement cycles test requirement per EN1060-3, part 8.10. This warranty extends only to consumers and does not extend to Retailers. To obtain warranty service on your Homedics product, contact a Consumer Relations Representative. Please make sure to have the model number of the product available. Homedics does not authorize anyone, including, but not limited to, Retailers, the subsequent consumer purchaser of the product from a Retailer or remote purchasers, to obligate Homedics in any way beyond the terms set forth herein. This warranty does not cover damage caused by misuse or abuse, accident, the attachment of any unauthorized accessory; alteration to the product; improper installation; unauthorized repairs or modifications; improper use of electrical/power supply; loss of power; dropped product; malfunction or damage of an operating part from failure to provide manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of Homedics. This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adoption to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. THERE SHALL BE NO OTHER WARRANTIES EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS OR ANY OTHER OBLIGATION ON THE PART OF THE COMPANY WITH RESPECT TO PRODUCTS COVERED BY THIS WARRANTY. HOMEDICS SHALL HAVE NO LIABILITY FOR ANY INCIDENTAL, CONSEQUENTIAL OR SPECIAL DAMAGES. IN NO EVENT SHALL THIS WARRANTY REQUIRE MORE THAN THE REPAIR OR REPLACEMENT OF ANY PART OR PARTS WHICH ARE FOUND TO BE DEFECTIVE. WITHIN THE EFFECTIVE PERIOD OF THE WARRANTY, NO REFUNDS WILL BE GIVEN. IF REPLACEMENT PARTS FOR DEFECTIVE MATERIALS ARE NOT AVAILABLE, HOMEDICS RESERVES THE RIGHT TO MAKE PRODUCT SUBSTITUTIONS IN LIEU OF REPAIR OR REPLACEMENT.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resaled products, including but not limited to sale of such products on internet auction sites and/or sales of such products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate as to any products or parts thereof which are repaired, replaced, altered, or modified, without the prior express and written consent of Homedics. This warranty provides you with specific legal rights. You may have additional rights which may vary from state to state. Because of individual state regulations, some of the above limitations and exclusions may not apply to you.

For more information regarding our product line in the USA, please visit: [www.homedics.com](http://www.homedics.com)

©2023 Homedics, LLC. All rights reserved.  
Homedics and the Smart Measure word mark and logo are registered trademarks of Homedics, LLC. American Heart Association is a registered trademark of the American Heart Association, Inc. The Bluetooth® word mark and logo are registered trademarks owned by Bluetooth SIG, Inc and any use of such marks by Homedics, LLC is under license. The App Store is a service mark of Apple Inc. Google Play and the Google Play logo are trademarks of Google LLC.  
IB-BPW800BTWTB  
P/N: 323103709 VER:002

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DISTRIBUIDO por: Homedics, LLC