## PULSE OXIMETER

### INSTRUCTIONS FOR USE

Model: CMS50D1A  
P/N: PX-131COB  
Doc# L-03969, Rev. 1

#### INTRODUCTION

In case of modifications and software upgrades, the information contained in this document is subject to change without notice. The Manual describes, in accordance with the Pulse Oximeter's features and requirements, main structure, functions, specifications, correct methods for transportation, installation, usage, operation, repair, maintenance and storage, etc., as well as the safety procedures to protect both the user and equipment. Refer to the respective sections for details.

Please read the Manual very carefully before using this equipment. These instructions describe the operating procedures that must be followed. Failure to follow these instructions can cause measuring abnormality, equipment damage and personal injury.

The manufacturer is NOT responsible for the safety, reliability and performance issues and any monitoring abnormality, personal injury and equipment damage due to user's negligence of the operation instructions. The manufacturer's warranty service does not cover such faults. Due to product improvements, the specific products you received may have slight differences to the description in this User Manual.

This product can be used repeatedly, with a life expectancy of 3 years.

#### WARNING

An uncomfortable or painful feeling may appear if using the device ceaselessly, especially for the microcirculation barrier users. It is recommended that the sensor should not be applied to the same finger for over 2 hours.

- Use the pulse oximeter only for its intended use. DO NOT use on a finger that has an edema and tender tissue.
- The light (the infrared is invisible) emitted from the device is harmful to the eyes. DO NOT stare into the light.
- User cannot use enamel or other makeup on the finger when using.
- User's fingernail cannot be too long.
- Please read the relative content about the clinical restrictions and cautions.
- This device is not intended for treatment or diagnosis.

#### IMPORTANT PRODUCT NOTICES AND SAFETY INSTRUCTIONS

- Check the main unit and all accessories periodically to make sure that there is no visible damage that may affect user's safety and monitoring performance about cables and transducers. It is recommended that the device should be inspected once a week at least. When there is obvious damage, stop using the oximeter.
- DO NOT disassemble this device. If in need of repair, refer to the warranty section of this manual.
- The oximeter cannot be used together with devices not specified in the User's Manual.
- Use only the accessories recommended in this manual, or included in the box.
- Under normal usage conditions, this device does not require calibration or technical maintenance other than infrequent replacement of the battery.

#### WARNINGS:

- Explosive hazard - DO NOT use the oximeter in an environment with inflammable gas such as some ignitable anesthetic agents.
- A person with allergies to rubber cannot use this device.
- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- When using your pulse oximeter, basic precautions should always be followed. Please read and follow all instructions and warnings before using this product. Save these instructions for future reference.
- Please don't measure this device with function test paper for the device's related information.
- DO NOT service or maintain this device while in use.
- The user is an intended operator.
- The probe of the device is the applied part.

#### ATTENTION:

- Keep the oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.
- If the oximeter gets wet, please stop operating it.
- When it is carried from cold environment to warm or humid environment, please do no use it immediately.
- DO NOT operate keys on front panel with sharp materials.
- High temperature or high pressure steam disinfection of the oximeter is not permitted.
- Refer to the "Care and Maintenance" section for instructions on cleaning and disinfection.
- DO NOT immerse the unit in water as this will result in damage to the unit.
- When cleaning the device with water, the temperature should be lower than 140°F (60°C).
- As to the fingers which are too thin or too cold, it would probably affect the normal measure of the users' SpO<sub>2</sub> and pulse rate, please clip the thick finger such as thumb and middle finger deeply enough into the probe.
- DO NOT use the device on infant or neonatal users.
- The product is suitable for adults weighing between 88 lbs to 242.5 lbs (40 kg to 110 kg).
- The device may not work for all users. If you are unable to achieve stable readings, discontinue use.
- The update period of data is less than 5 seconds, which will change according to the user's pulse rate.
- The waveform is normalized. Please read the measured value when the waveform on screen is equally and steady-going, here this measured value is optimal value, and the waveform at the moment is the standard one.
- If some abnormal conditions appear on the screen while using, pull out the finger and reinsert to restore normal use.
- The lanyard attached to the product is made from non-allergic material, if particular groups are sensitive to the lanyard, stop using it. DO NOT wear around the neck to avoid risk of injury to the neck.
- The instrument does not have a low-voltage alarm function, it only shows the low-voltage, please change the battery when the battery energy is used up.
- The instrument does not have an alarm function. Do not use the device in situations where alarms are required.
- Batteries must be removed if the device is going to be stored for more than one month, or else the batteries may leak.
- A flexible circuit connects the two parts of the device. DO NOT twist or pull on the connection.

#### INDICATION FOR USE:

The pulse oximeter is a non-invasive device intended for the spot-check of oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and pulse rate of adults through the finger in home use environments. For use with sporting and aviation activities, and to monitor heart rate during exercise. Not intended for single use and out-of-hospital transport use.

#### OVERVIEW

The pulse oxygen saturation is the percentage of HbO<sub>2</sub> in the total Hb in the blood, so-called the O<sub>2</sub> concentration in the blood. It is an important bio-parameter for the respiration. For the purpose of measuring the SpO<sub>2</sub> more easily and accurately, our company developed the Pulse Oximeter. At the same time, the device can measure the pulse rate simultaneously.

The Pulse Oximeter features in small volume, low power consumption, convenient operation and being portable. It is only necessary for user to put one of his fingers into a fingertip photoelectric sensor for measurement, and the display screen will directly show measured value of Hemoglobin Saturation.

#### FEATURES:

- Operation of the product is simple and convenient.
- The product is small in volume, light in weight and convenient to carry.
- Power consumption of the product is low and the 2 originally equipped AAA batteries can be operated continuously for 20 hours.
- The product will enter standby mode when no signal is in the product within 5 seconds.
- Display direction changes automatically for easy viewing.

#### MAJOR APPLICATIONS AND SCOPE OF APPLICATION

The Pulse Oximeter can be used in measuring pulse oxygen saturation and pulse rate through the finger. The product is suitable for family use (It can be used before or after doing sports, and it is not recommended to use the device during the process of doing sports).

The problem of overrating would emerge when the user is suffering from toxicosis which is caused by carbon monoxide, the device is not recommended to be used under this circumstance.

#### ENVIRONMENT REQUIREMENTS

##### Storage Environment

- Temperature: -40°F (-40°C) ~ +140°F (+60°C)
- Relative humidity: ≤95%
- Atmospheric pressure: 500 hPa ~ 1060 hPa

##### Operating Environment

- Temperature: 50°F (10°C) ~ 104°F (40°C)
- Relative Humidity: ≤75%
- Atmospheric pressure: 700 hPa ~ 1060 hPa

#### PRINCIPLE AND CAUTION

##### Principle of Measurement

Principle of the Oximeter is as follows: An experience formula of data process is established taking use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive Hemoglobin (Hb) and Oxyhemoglobin (HbO<sub>2</sub>) in glow & near-infrared zones. Operation principle of the instrument is: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning & Recording Technology, so that two beams of different wavelength of lights can be focused onto human nail tip through perspective clamp finger-type sensor. Then measured signal can be obtained by a photosensitive element, information acquired through which will be shown on screen through treatment in electronic circuits and microprocessor.


Figure 1. Operating Principle

##### Caution

1. The finger should be placed properly (refer to Figure 7), or else it may cause inaccurate measurement.
2. The SpO<sub>2</sub> sensor and photoelectric receiving tube should be arranged in a way with the subject's arteriole in a position there between.
3. The SpO<sub>2</sub> sensor should not be used at a location or limb tied with arterial canal or blood pressure cuff or receiving intravenous injection.
4. Make sure the optical path is free from any optical obstacles like rubberized fabric.
5. Excessive ambient light may affect the measuring result. It includes fluorescent lamp, dual ruby light, infrared heater, direct sunlight, etc.
6. Strenuous action of the subject or extreme electrosurgical interference may also affect the accuracy.
7. User cannot use enamel or other makeup on the finger when using.

#### CLINICAL RESTRICTIONS

1. As the measurement is taken on the basis of arteriole pulse, substantial pulsating blood flow of subject is required. For a subject with weak pulse due to shock, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO<sub>2</sub> waveform (PLETH) will decrease. In this case, the measurement will be more sensitive to interference.
2. For those with a substantial amount of staining dilution drug (such as methylene blue, indigo green and acid indigo blue), or carbon monoxide hemoglobin (CO-Hb), or methionine (Me+Hb) or thiosalicylic hemoglobin, and some with icterus problem, the SpO<sub>2</sub> determination by this monitor may be inaccurate.
3. The drugs like dopamine, procaine, prilocaine, lidocaine and butacaine may also be a major factor blamed for serious error of SpO<sub>2</sub> measure.
4. As the SpO<sub>2</sub> value serves as a reference value for judgement of anemic anoxia and toxic anoxia, some users with serious anemia may also report good SpO<sub>2</sub> measurement.

#### TECHNICAL SPECIFICATIONS

1. Display Format: LCD Display;  
SpO<sub>2</sub> Measuring Range: 0% ~ 100%;  
Pulse Rate Measuring Range: 30 bpm ~ 250 bpm;  
Pulse Wave Display: bar-graph display and the waveform display.
2. Power Requirements: 2 x 1.5 V AAA alkaline battery (or use rechargeable batteries), adaptable range: 2.6 V ~ 3.6 V.
3. Power Consumption: Less than 50 mA.
4. Resolution: 1% for SpO<sub>2</sub> and 1 bpm for Pulse Rate.
5. Measurement Accuracy: ±2% in stage of 70% ~ 100% SpO<sub>2</sub>, and meaningless when stage being smaller than 70%. ±2 bpm or ±2% (select larger) for Pulse Rate. Clinical Trial: SpO<sub>2</sub> regression plot and Bland-Altman plot, Refer to Figure 2 & Figure 3.
6. Measurement Performance in Weak Filling Condition: SpO<sub>2</sub> and pulse rate can be shown correctly when pulse-filling ratio is 0.4%. SpO<sub>2</sub> error is ±4%, pulse rate error is ±2 bpm or ±2% (select larger).
7. Resistance to surrounding light: The deviation between the value measured in the condition of man-made light, indoor natural light and darkroom is less than ±1 %.
8. It is equipped with a switch function: The product will enter standby mode when no signal is in the product within 5 seconds.
9. Optical Sensor: Red light (wavelength is 660 nm, 6.65 mW), infrared (wavelength is 905 nm, 6.75 mW).


Figure 2. SpO<sub>2</sub> regression plot


Figure 3. Bland-Altman plot

#### PRODUCT SET UP


Figure 4. Front View


Figure 5. Attaching the lanyard

#### ATTACHING THE LANYARD

1. Put the end of the lanyard through the hole.
2. Put another end of the lanyard through the first one and then tighten it.

#### BATTERY INSTALLATION

1. Refer to Figure 6. and insert the 2 AAA size batteries properly in the right direction.
2. Put back the cover.

**⚠ Please take care when you insert the batteries for the improper insertion may damage the device.**

- Replace all batteries at one time (as simultaneous set). Use only 1.5V AAA alkaline batteries. Use only the size and type of batteries specified.
- When installing batteries, observe proper +/- polarities. Incorrect installation of batteries may cause damage to the unit.
- Clean contacts on battery and in battery compartment with a soft, dry cloth each time you install batteries.
- Batteries are hazardous waste. DO NOT dispose of them together with household garbage.
- DO NOT dispose of batteries in fire. Batteries may explode or leak.
- Recycle or dispose of properly in accordance with local, state, province, and country regulations.
- Keep batteries away from children.
- Batteries are small objects and could be ingested. If swallowed, contact a physician at once.
- Remove batteries promptly from the device when the batteries are depleted.


Figure 6. Battery Installation

## OPERATING GUIDE

Figure 7. Finger Position

- With the batteries in the device press the back of the device to open the clip.
- Place the finger between the rubber cushions of the clip as shown in Figure 7 (make sure the finger is in the right position), and then clip the finger.
- Press the power button to start measurement.
- DO NOT move the finger or body during measurement.
- The measurement information will be shown directly on the display.
- Long press the button to adjust brightness preference. Choose from 5 levels of brightness. When in Standby mode, a short press of the button will exit it.
- Display viewing angle will change automatically according to the direction of the hand.

**⚠ Fingernails and the luminescent tube should be on the same side.**

## CARE AND MAINTENANCE

- Replace the batteries if the low battery symbol appears on the display, or when any button is pressed and nothing is displayed on the screen.
- Please clean the surface of the device before using. Wipe the device with rubbing alcohol first, let it air dry, or clean it using a soft, dry piece of clean fabric.
- Use rubbing alcohol to disinfect the product after use to prevent spreading infection between uses.
- Remove batteries when unit is not in use for extended periods of time.
- The packed device can be transported by ordinary conveyance or according to transport contract. The device cannot be transported mixed with toxic, harmful, corrosive material.
- Storage environment of the device is -40°F (-40°C) to +140°F (60°C), < 95% RH.

**⚠ High-pressure sterilization cannot be used on the device. DO NOT immerse the device in liquid.**  
It is recommended that the device should be kept in a dry environment.  
Humidity may reduce the useful life of the device, or even damage it.

## TROUBLESHOOTING

PROBLEM	CAUSE	SOLUTION
The SpO <sub>2</sub> and pulse rate are not displayed normally	1. The finger is not properly positioned. 2. The user's SpO <sub>2</sub> is too low to be detected.	1. Place the finger properly and try again. 2. Try again. Go to a hospital for a diagnosis if you are sure the device works all right.

The SpO <sub>2</sub> and pulse rate are not displayed normally	1. The finger is not placed inside deep enough. 2. The finger is shaking or the user is moving.	1. Place the finger properly and try again. 2. Keep the user calm.
The device cannot be turned on	1. The batteries are drained or almost drained. 2. The batteries are not inserted properly. 3. The device is malfunctioning.	1. Change batteries. 2. Reinstall batteries. 3. Please contact Consumer Relations.
The display is off suddenly	1. The product will enter standby mode when no signal is detected within 5 seconds. 2. The batteries are almost drained.	1. Normal. 2. Change batteries.

## SYMBOLS

SYMBOL	DESCRIPTION	SYMBOL	DESCRIPTION
	Type BF		WEEE (2002/96/EC)
	Refer to instruction manual/booklet		Ingress of liquids rank
%SpO <sub>2</sub>	The pulse oxygen saturation(%)		Manufacturer
PR bpm	Pulse rate (bpm)		Manufacture Date
	Full-voltage		Storage and Transport Temperature limitation -40°C ~ +60°C
	The battery voltage indication is deficient (change the battery in time avoiding the inexact measure)		Storage and Transport Humidity limitation 5% ~ 95%
	1. No finger inserted 2. An indicator of signal inadequacy		Storage and Transport Atmospheric pressure limitation
	Battery positive electrode		This side UP
	Battery negative electrode		Fragile, handle with care
	1. Exit standby mode 2. Change brightness of the screen		Keep dry
SN	Serial number		Recyclable
	Alarm inhibit		

## SPECIFICATIONS

DISPLAY INFORMATION	DISPLAY MODE
The Pulse Oxygen Saturation(SpO <sub>2</sub> )	LCD
Pulse Rate(PR)	LCD
Pulse Intensity (bar-graph)	LCD bar graph display
Pulse wave	LCD
SpO <sub>2</sub> Parameter Specification	
Measuring range	0% ~ 100% (the resolution is 1%).
Accuracy	70% ~ 100%: ±2%, below 70% unspecified.
Optical Sensor	Red light (wavelength is 660 nm) Infrared (wavelength is 905 nm)
Pulse Parameter Specification	
Measuring range	30 bpm ~ 250 bpm (the resolution is 1 bpm).
Accuracy	±2 bpm or ±2% select larger
Pulse Intensity	
Range	Continuous bar-graph display, the higher display indicate the stronger pulse.
Battery Requirement	
1.5V (AAA size) alkaline batteries x 2 or rechargeable battery	
Battery Life	
Up to 20 hours	
Dimensions and Weight	
Dimensions	2.36"(L) X 1.42"(W) X 1.26"(H) 60(L) X 36(W) X 32(H) mm
Weight	About 0.002 lbs (57g) with the batteries

## APPENDIX

Guidance and manufacturer's declaration-electromagnetic emission for all EQUIPMENT and SYSTEMS
<b>GUIDANCE AND MANUFACTURE'S DECLARATION—ELECTROMAGNETIC EMISSION</b>

The PX-131CO Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the PX-131CO Pulse Oximeter should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR II	Group 1	The PX-131CO Pulse Oximeter uses RF energy only for their internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR II	Class B	The PX-131CO Pulse Oximeter is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emission IEC 61000-3-3	Not applicable	

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM for EQUIPMENT or SYSTEM that not LIFE-SUPPORTING

### RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE PX-131CO PULSE OXIMETER

The PX-131CO Pulse Oximeter is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the PX-131CO Pulse Oximeter can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the PX-131CO Pulse Oximeter as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
d = $\left[ \frac{3.5}{E_1} \right] \sqrt{P}$	d = $\left[ \frac{3.5}{E_1} \right] \sqrt{P}$	d = $\left[ \frac{7}{E_1} \right] \sqrt{P}$	d = $\left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## 2-YEAR LIMITED WARRANTY

Homedics sells its products with the intent that they are free of defects in manufacture and workmanship for a period of 2 years from the date of original purchase, except as noted below. Homedics warrants that its products will be free of defects in material and workmanship under normal use and service. This warranty extends only to Consumers and does not extend to Retailers. To obtain warranty service on your Homedics product, contact a consumer relations representative for assistance. Please make sure to have the model number of the product available.

Homedics does not authorize anyone, including, but not limited to, Retailers, the subsequent consumer purchaser of the product from a Retailer or remote purchasers, to obligate Homedics in any way beyond the terms set forth herein. This warranty does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; improper installation; unauthorized repairs or modifications; improper use of electrical power/supply; loss of power; dropped product; malfunction or damage of an operating part from failure to provide manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of Homedics.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adoption to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized or repair of products damaged by these modifications is not covered under this warranty.

**THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. THERE SHALL BE NO OTHER WARRANTIES EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL THE PARTIES INVOLVED IN THE PURCHASE OR SALE OF THIS PRODUCT BE LIABLE FOR INCIDENTAL OR SPECIAL DAMAGES. IN NO EVENT SHALL THIS WARRANTY REQUIRE MORE THAN THE REPAIR OR REPLACEMENT OF ANY PART OR PARTS WHICH ARE FOUND TO BE DEFECTIVE WITHIN THE EFFECTIVE PERIOD OF THE WARRANTY. NO REFUNDS WILL BE GIVEN IF REPLACEMENT PARTS FOR DEFECTIVE MATERIALS ARE NOT AVAILABLE. HOMEDICS RESERVES THE RIGHT TO MAKE PRODUCT SUBSTITUTIONS IN LIEU OF REPAIR OR REPLACEMENT.**

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resold products, including but not limited to the sale of such products on internet auction sites and/or sales of such products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate as to any products or parts thereof which are repaired, replaced, altered, or modified, without the prior express and written consent of Homedics.

This warranty provides you with specific legal rights. You may have additional rights which may vary from state to state. Because of individual state regulations, some of the above limitations and exclusions may not apply to you.

For more information regarding our product line in the USA, please visit: [www.homedics.com](http://www.homedics.com).

For service or repair, do not return this unit to the retailer. Contact Homedics Consumer Relations:

Email: [cservice@homedics.com](mailto:cservice@homedics.com)  
Phone: 1-800-466-3342

Business Hours: 8:30am-7pm EST, Monday-Friday

Homedics® is a registered trademark of Homedics, LLC.

©2023 Homedics, LLC. All rights reserved.

IB-PX131COC

X.XXXXXXX

In USA Distributed by:  
Homedics USA, LLC  
3000 N Pontiac Trail  
Commerce Township, MI 48390

Printed in China

NOTE 1: At 80MHz and 800MHz, the higher frequency range applies.  
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the PX-131CO Pulse Oximeter is used exceeds the applicable RF compliance level above, the PX-131CO Pulse Oximeter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the PX-131CO Pulse Oximeter.

b. Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3V/m.

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Homedics® es una marca registrada de Homedics, LLC.